

Ferno iNX Stretcher

Securing the Patient

1. Unbuckle all straps in preparation for the patient.
2. Position patient in center of stretcher with patient's back flat, and buttocks at the joint of the stretcher.
3. Lay shoulder straps on the patient's chest with the red webbing and metal links extended to the patient's waist.
4. Take the pelvis strap buckle (male part) and guide it through the slots in the shoulder strap links. Buckle the pelvis strap to connect all red parts.
5. Adjust the shoulder and pelvis straps to be snug on the patient.
6. Fasten the black chest strap across the patient's rib cage.
7. Fasten the black leg restraint over the patient's legs and snug.

3.



4.



6.



* Some older Ferno stretchers may use a different shoulder restraint. Refer to the appropriate Ferno manual and videos for complete instructions and information on product storage and maintenance.

Stryker Power Pro XT Stretcher (With XPR restraints*)

Securing the Patient

1. Unbuckle all straps in preparation for the patient.
2. Position patient in center of stretcher with patient's back flat, and buttocks at the joint of the stretcher.
3. Lay one shoulder strap diagonally across the patient's chest, buckle, and snug the harness. Repeat with other shoulder strap crossing the chest.
4. Position the waist restraint over the patient's pelvic bones. Buckle it and snug the strap.
5. Place the leg restraint over the upper legs. Buckle and snug the straps.
6. Fasten the second leg restraint over the patient's lower legs. Buckle and snug the strap.



3.



4.



5. & 6.



* Some older Stryker stretchers may use a different shoulder restraint. Refer to the appropriate Stryker manual and videos for complete instructions and information on product storage and maintenance.